

Holds Six 2.5-Liter Bottles, 2 Doors, Manual Close, Wood Laminate Corrosives Safety Cabinet, Countertop, Blue - 24120

Product Images



Short Description

- Durable and highly chemical-resistant
- Safely store highly corrosive liquids
- Ideal for smaller quantities
- Long-lasting wood construction
- Haz-Alert™ reflective labels

Description

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous

explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.

Safely store highly corrosive liquids such as phenol, nitric acid and sulfuric acid in these wood laminate storage cabinets. Ideal for smaller quantities used at workbench level. Long lasting wood construction along with stainless steel handles and hinges provide exceptional resistance to corrosion from harsh acids. Shelf adjusts on 1-1/4" (32-mm) centers to accommodate a wide variety of containers. Cabinets include adjustable leg levelers, Haz-Alert™ reflective labels, and double key set for safety and security.

Specifications

| | |
|-------------------------|-----------------------------|
| Model No | 24120 |
| UPC | 697841006607 |
| International Model No. | 24120 |
| Color | Blue |
| Material Specifications | Wood laminate |
| Capacity | Six 2½-Liter Bottles |
| Liter Capacity | Six 2½-Liter Bottles |
| Number of Doors | 2 |
| Adjustable Shelves | Yes |
| Dimensions, Exterior | 18.5" H x 24" W x 16" D |
| | 470mm H x 610mm W x 406mm D |
| Dimensions, Interior | 16" H x 22.5" W x 15.2" D |
| | 406mm H x 571mm W x 387mm D |
| Height (Inches) | 18.5 |
| Width (Inches) | 24 |
| Depth | 16 |
| Net Weight, lbs | 45.0 |
| Net Weight, kg | 20.4 |