

Holds 49, 2.5-Liter Bottles, 2 Shelves, 2 Doors, Manual Close, Wood Laminate Corrosives Safety Cabinet, Blue - 24150

Product Images



Short Description

- Durable and highly chemical-resistant
- Safely store highly corrosive liquids
- Long-lasting wood construction
- Haz-Alert™ reflective labels

Description

Research environments using many different types of chemicals have special storage requirements. Whether using acids, bases or solvents, it is critical that incompatible liquids are segregated to avoid disastrous explosions or fire. The SDS (Safety Data Sheet) should be reviewed to determine chemical characteristics and recommended storage practices.

Safely store highly corrosive liquids such as phenol, nitric acid and sulfuric acid in these wood laminate storage cabinets. Long lasting wood construction along with stainless steel handles and hinges provide exceptional resistance to corrosion from harsh acids. Shelf adjusts on 1-1/4" (32-mm) centers to accommodate a wide variety of containers. Cabinets include adjustable leg levelers, Haz-Alert™ reflective labels, and double key set for safety and security.

Specifications

Model No	24150
UPC	697841006638
International Model No.	24150
Color	Blue
Pallet Quantity	2
Material Specifications	Wood laminate
Liter Capacity	Forty-nineThirty six 2½-Liter Bottles
Number of Doors	2
Number of Shelves	2
Adjustable Shelves	Yes
Dimensions, Exterior	60" H x 42" W x 17.875" D 1524mm H x 1067mm W x 454mm D
Dimensions, Interior	54.5" H x 40.5" W x 17.1" D 1385mm H x 1028mm W x 435mm D
Height (Inches)	60
Width (Inches)	42
Depth	17.875
Net Weight, lbs	207.0
Net Weight, kg	93.9