

# 1 Gallon, Oval Polyethylene Safety Shelf Can, Stainless Steel Hardware, Bottom Self-Close S/S Faucet, Red - 14169

## Product Images



## Short Description

- Non-metallic laboratory safety can
- Self-closing faucet
- Features pressure-relief cap with stainless steel flame arrester
- FM approved

## Description

---

Nonmetallic Laboratory Safety Cans have all the benefits of Justrite's Type I polyethylene safety cans, with the addition of a self-closing faucet for easy dispensing. All fill spouts have a pressure-relief cap with stainless steel flame arrester. Hardware featured on the can top is stainless steel for maximum corrosion resistance and long life. The 5-gallon size also features a vacuum-relief valve in the lid allowing liquids to flow smoothly. Justrite's unique patented, current-carrying carbon insert embedded into the container provides for proper grounding.

**Nonmetallic Shelf Cans** have a faucet mounted 2 inches above the bottom to eliminate lifting the cans during dispensing. FM approved.

## Specifications

---

Model No	14169
UPC	697841005709
International Model No.	14169
FM Approval	Yes
OSHA Compliance	Yes
Complies with NFPA Code 30	Yes
Color	Red
Material Specifications	Polyethylene
Dimensions, Exterior	11.5" H x 7.25" Outer Dia. 292mm H x 184mm Outer Dia.
Capacity	1 gallon
Liter Capacity	4
Fill/Pour Hole inch	1.45
Fill/Pour Hole mm	37
Net Weight, lbs	2.2
Net Weight, kg	1.0