

Coalescing/Carbon Filter With Disconnect, for HPLC Disposal Cans With Stainless Steel Disconnect - 28162

Product Images



Short Description

- Minimize odors
- Easy to install
- Meet EPA 40 CFR 264.173

Description

Coalescing/Carbon Filters minimize odor by capturing 99% of liquid particulates released in the form of gaseous vapor. Vent filter is easy to install, snapping conveniently into the “Quick-Disconnect” fitting. A color-changing Carbon Cartridge (top portion) changes from purple to black, as hydrocarbons are absorbed indicating replacement is necessary. The Coalescing Filter (bottom portion) can be reused several times (three or four carbon filter cycles) before replacing. Filters have been designed for use with all two and five-gallon TF, PP and custom-built “Quick-Disconnect” cans.

Meet EPA 40 CFR 264.173.

Specifications

Model No	28162
UPC	697841074644
International Model No.	28162
EPA Compliance	Yes
Color	white
Material Specifications	Plastic
	89mm Outer Dia.
Net Weight, lbs	1.0
Net Weight, kg	0.5