

Replacement 1/16" OD Tube Fittings With PCTFE Ferrules, for HPLC Stainless Steel Manifold, Set of 4 - 28174

Product Images



Short Description

- HPLC Stainless Steel Manifold Kit, Tube Fittings with PCTFE Ferrules

Description

HPLC Stainless Steel Manifold fittings, four each: 1/16" (2-mm) OD tube fittings with corresponding PCTFE ferrules for compression fit of waste line.

Specifications

Model No	28174
UPC	697841097223
International Model No.	28174
Material Specifications	Polypropylene
Dimensions, Exterior	0.063" Outer Dia.
Net Weight, lbs	1.5
Net Weight, kg	0.7