

1/16" OD Hard-Wall Tubing Compression Fittings With Ferrules for HPLC Poly Manifold, Package of 6 - 28188

Product Images



Short Description

- Sold in sets of six
- Include ferrules for compression fit to waste line

Description

The compression fittings, sold in sets of six, accommodate hard-wall tubing size of 1/16-inch OD. The fittings are approximately 1" (25-mm) long and made of ETFE (ethylene tetrafluoroethylene) — a fluorine based plastic. They include ferrules for a compression fit to the waste line.

Specifications

Model No	28188
UPC	697841121461
International Model No.	28188
Color	Blue
Material Specifications	Resin
Dimensions, Exterior	0.75" L x 0.25" Outer Dia.
Net Weight, lbs	0.1
Net Weight, kg	0.0