


5 Gallon Capacity, HPLC Can Spill Basin, Polyethylene, Translucent - 84003

Product Images



Short Description

- Keeps work areas clean around HPLC cans
- Captures leakage from HPLC inflow tubs
- Protects cans stored on floor against accidental bumping

Description

The durable, polyethylene Spill Basin keeps work areas clean around the HPLC cans, and captures any leakage from HPLC inflow tubes. The basin provides additional protection from accidental “bumping” when cans are stored on the floor. Sized to accept Justrite’s two and five-gallon poly disposal cans with a 12”

(305mm) diameter, this leakproof basin makes cleanup easy and offers a spill capacity of up to five gallons in accordance with EPA requirements. The translucent color detects spills even at a distance.

Specifications

Model No	84003
UPC	697841124332
International Model No.	84003
Pallet Quantity	48
Color	white
Material Specifications	Polyethylene
Dimensions, Exterior	7.625" H x 14.875" W x 14.875" D
	194mm H x 378mm W x 378mm D
Capacity	5 gallon
Liter Capacity	19
Net Weight, lbs	1.5
Net Weight, kg	0.7