


## 5 Gallon Capacity, HPLC Can Spill Basin, Polyethylene, Translucent - 84003

### Product Images

---



### Short Description

---

- Keeps work areas clean around HPLC cans
- Captures leakage from HPLC inflow tubs
- Protects cans stored on floor against accidental bumping

### Description

---

The durable, polyethylene Spill Basin keeps work areas clean around the HPLC cans, and captures any leakage from HPLC inflow tubes. The basin provides additional protection from accidental “bumping” when cans are stored on the floor. Sized to accept Justrite’s two and five-gallon poly disposal cans with a 12”

(305mm) diameter, this leakproof basin makes cleanup easy and offers a spill capacity of up to five gallons in accordance with EPA requirements. The translucent color detects spills even at a distance.

## Specifications

---

Model No	84003
UPC	697841124332
Brand	Justrite
International Model No.	84003
Pallet Quantity	48
Color	White
Material Specifications	Polyethylene
Dimensions, Exterior	7.625" H x 14.875" W x 14.875" D
Dimensions, Exterior	194mm H x 378mm W x 378mm D
Capacity	5 gallon
Liter Capacity	19
Net Weight, lbs	1.5
Net Weight, kg	0.7