


44.375" W x 18" D Steel Shelf for 31 Gallon Under Fume Hood Safety Cabinet, ChemCor®, SpillSlope® - 29948

Product Images



Short Description

- Customize ChemCor® cabinets with additional shelves
- Provide superior protection against corrosion and off-gassing
- Patented SpillSlope® galvanized steel shelves direct spills to back and bottom of a leakproof sump

- Fits 31-gallon Under Fume Hood Cabinet

Description

Customize your ChemCor® cabinet with additional shelves

Shelves have the ChemCor® thermoplastic coating and provide superior protection against corrosion from drips and off-gassing from harsh chemicals

Patented SpillSlope® galvanized steel shelves direct spills to back and bottom of a leakproof sump.

Fits 31-gallon Under Fume Hood Cabinet.

Specifications

Model No	29948
Accessory Type	Shelves and Trays
UPC	697841130265
Brand	Justrite
International Model No.	29948
Color	Gray
Material Specifications	Steel
Load Capacity Compliance	350 lb. (159 kg.) uniformly distributed load (UDL)
Adjustable Shelves	Yes
Dimensions, Exterior	44.375" W x 18" D x 1.428" H
Dimensions, Exterior	1127mm W x 457mm D
Net Weight, lbs	13.0
Net Weight, kg	5.9