

## Filter, Carbon, for VaporTrap Solvent Waste Kits, Package of 1 - #12848

### Product Images

---



### Short Description

---

- Prevent harmful HPLC vapors from being released into the air
- Activated Carbon filter prevents volatile organic compound vapors from being released into the air
- Protect workers from inhaling dangerous fumes

### Description

---

Designed to be used in conjunction with VaporTrap™ solvent waste systems, these carbon exhaust filters trap potentially harmful vapors emitted from solvents in HPLC applications, protecting workers from inhaling dangerous fumes.

Depending on application, filter life may vary. Optional filter replacement indicator makes it easy to see when

filter is saturated and requires replacement.

## Specifications

---

Model No	12848
UPC	697841181915
International Model No.	12848
Material Specifications	Charcoal Activated Carbon
Dimensions, Exterior	7" L x 1.9" Outer Dia. 178mm L x 48mm Outer Dia.
Net Weight, lbs	0.5
Net Weight, kg	0.2