

## Fittings for Carboy Adapters, six 1/8" fittings and six 1/4" plugs, Package of 1 - #12851

### Product Images

---



### Short Description

---

- Solvent waste kit replacement fittings

### Description

---

Replacement fittings for VaporTrap™ solvent waste system adapter includes six 1/8-inch nut and ferrule fittings, six 1/4-inch fitting plugs, one straight filter fitting, one 90-degree elbow fitting, and one filter fitting plug.

## Specifications

---

Model No	12851
UPC	697841181946
International Model No.	12851
Dimensions, Exterior	1" H x 1" W x 1" L
	25mm H x 25mm W x 25mm L
Cap Size	70mm
Net Weight, lbs	0.2
Net Weight, kg	0.1