

53mm Carboy Cap, Open Top, Closed Adapter - 12891

Product Images



Short Description

- Aseptic handling solutions for high purity fluids and sensitive biologicals
- Eliminate open-air pouring of high valued products with this leak proof, closed system

- High purity platinum-cured silicone insert with or without tubing

Description

- Eliminate open-air pouring of high valued products with this leak proof, closed system
- High purity platinum-cured silicone insert with or without tubing

Designed for aseptic fluid transfer, sampling, or sealing of high purity fluids and sensitive biologicals in pharmaceutical and laboratory applications, these systems features a platinum-cured silicone insert, with or without silicone tubing, and a polypropylene cap that fits bottles or carboys with 53-mm, 83-mm, or 120-mm closures. All components can be autoclaved at 15 PSI and 250°F (121°C) for 20 minutes.

Note: Height includes pre-mounted tubing.

Specifications

| | |
|-------------------------|--|
| Model No | 12891 |
| UPC | 697841181168 |
| International Model No. | 12891 |
| Material Specifications | Polypropylene |
| Dimensions, Exterior | 1.8" H x 3" Outer Dia. 46mm H x 76mm Outer Dia. |
| Net Weight, lbs | 0.8 |
| Net Weight, kg | 0.4 |