
13.5 Liter HDPE Carboy, UN/DOT Certified, 70mm Cap - 12950

Product Images



Short Description

- UN/DOT approved carboy ensures safe transport of hazardous liquids
- Meets UN 3H1 standards
- High-density polyethylene (HDPE) construction for excellent chemical compatibility

Description

- UN/DOT approved carboy for packaging group II, and hazardous chemicals
- Meets UN 3H1 standards
- High-density polyethylene (HDPE) construction for excellent chemical compatibility

Designed for storage and secure shipments of liquid waste, these UN/DOT certified carboy (UN Rating

3H1/Y1.8/100) are designed for single-use and easy disposal of laboratory waste. The leakproof 70-mm HDPE cap includes a chemical resistant EPDM gasket. The ergonomic HDPE carboy features a large top handle for a full hand grip, and provides excellent chemical compatibility against aggressive solvents.

Specifications

Model No	12950
UPC	697841182110
International Model No.	12950
D.O.T Approval	Yes
UN Certification No. for Salvage drum	3H1/Y1.8/100
Material Specifications	High Density Polyethylene
Dimensions, Exterior	11" H x 9.25" W x 11" L 279mm H x 235mm W x 279mm L
Capacity	3.6 gallon
Liter Capacity	13.5
Net Weight, lbs	2.0
Net Weight, kg	0.9