
9 Gallon Spill Basin for 13.5 or 20 Liter Carboys, Plastic, Red - 12956

Product Images



Short Description

- Keep work area clean and capture leaks from HPLC tubes
- Capacity 9 gallons (34L)

Description

These secondary containers capture leaks and spills from the VaporTrap™ solvent waste systems and carboys. Constructed of a rugged high-density polyethylene, the spill basins are a perfect solution for keeping work areas clean and dry.

Specifications

Model No	12956
UPC	697841182158
International Model No.	12956
Material Specifications	High Density Polyethylene
Dimensions, Exterior	10" H x 15.5" W x 19.5" L
	254mm H x 394mm W x 495mm L
Capacity	9 gallon
Liter Capacity	34
Net Weight, lbs	3.1
Net Weight, kg	1.4