

# Type I Shelf Safety Can, 5 gallon, bottom faucet, Steel.

## Product Images

---



## Short Description

---

- Ideal for laboratory settings
- Self-close, pressure relief caps
- Vacuum-relief valve within the lid
- Faucet mounted 2 inches above the bottom
- FM approved

## Description

---

Justrite's Steel Laboratory Safety Cans have all the benefits of the Type I cans. Added features make these the choice of lab professionals for safely dispensing hazardous liquids into test tubes, small flasks and beakers. Fill spouts have self-close, pressure-relief caps and stainless steel flame arresters. These cans also feature a vacuum-relief valve within the lid for smooth liquid flow, and self-closing faucets for trouble-free dispensing. Cans with the smaller No. 08540 faucet offer easier dispensing in lab environments, while the No. 08902 faucet is better suited for heavy, industrial use.

**Shelf Cans** have a No. 08540 faucet mounted 2 inches above the bottom to eliminate lifting the cans during dispensing.

FM Approved

## Specifications

---

Model No	7150140 / 7150240
FM Approval	Yes
OSHA Compliance	Yes
Complies with NFPA Code 30	Yes
CARB Compliance	Yes
Material Specifications	Steel/Brass Faucet
Dimensions, Exterior	16.875" H x 11.75" Outer Dia. 429mm H x 298mm Outer Dia.
Capacity	5 gallon
Liter Capacity	19
Fill/Pour Hole inch	1.58
Fill/Pour Hole mm	40
Net Weight, lbs	6.0
Net Weight, kg	2.7

---

Color:

Red

Yellow

## Product Options