

Type I Shelf Safety Can, 5 gallon, bottom faucet, Steel.

Product Images



Short Description

- Ideal for laboratory settings
- Self-close, pressure relief caps
- Vacuum-relief valve within the lid
- Faucet mounted 2 inches above the bottom
- FM Approved

Description

Justrite's Steel Laboratory Safety Cans have all the benefits of the Type I cans. Added features make these the choice of lab professionals for safely dispensing hazardous liquids into test tubes, small flasks and beakers. Fill spouts have self-close, pressure-relief caps and stainless steel flame arresters. These cans also feature a vacuum-relief valve within the lid for smooth liquid flow, and self-closing faucets for trouble-free dispensing. Cans with the smaller No. 08540 faucet offer easier dispensing in lab environments, while the No. 08902 faucet is better suited for heavy, industrial use.

Shelf Cans have a No. 08902 faucet mounted 2 inches above the bottom to eliminate lifting the cans during dispensing.

FM Approved

Specifications

Model No	7150150 / 7150250
Brand	Justrite
FM Approval	Yes
OSHA Compliance	Yes
Complies with NFPA Code 30	Yes
Material Specifications	Steel/Brass Faucet
Dimensions, Exterior	16.875" H x 11.75" Outer Dia.
Dimensions, Exterior	429mm H x 298mm Outer Dia.
Capacity	5 gallon
Liter Capacity	19
Fill/Pour Hole inch	1.58
Fill/Pour Hole mm	40
Net Weight, lbs	8.0
Net Weight, kg	3.6

Color:

Red

Yellow

Product Options