


# Platinum-Cured Silicone O-Rings for Carboys, 120mm Caps, 5 Pack - 12961

## Product Images

---



## Short Description

---

- Replacement o-rings for carboy caps
- Fits 120mm caps
- Five to a package

## Description

---

Replacement o-rings for Carboy caps

## Specifications

---

Model No	12961
UPC	697841202757
Quantity Per Pack	5
Material Specifications	Platinum-Cured Silicone
Cap Size	Fits 120mm Carboy cap
Net Weight, lbs	1
Net Weight, kg	1