

Replacement 1/8" MNPT Black Fitting With 3/8" Hose Barb Plug, for VaporTrap™, 10-Pack - 12965

Product Images



Short Description

- Replacement Fittings, 1/8-in MNPT, 3/8-in Hose Barb Plug
- Convenient 10 pack

Description

The 1/8-in MNPT x 3/8" Hose Barb fittings pack is the replacement fittings for VaporTrap™ Solvent Waste System's fittings. It includes (10) 1/8-in MNPT x 3/8-in hose barb fittings. These tube fittings secure waste tubes to carboy providing the closed system necessary to prevent spills.

Specifications

Model No	12965
UPC	697841206236
Color	Black
Quantity Per Pack	10
Material Specifications	Polyethylene
MNPT	1/8-in
1/4" or 3/8" Tubing ID Hose Barb Ports	3/8-in
Net Weight, lbs	1
Net Weight, kg	1