

Replacement Frame for Lab Tilt Can 11273

Product Images



*Safety Can not included

Short Description

- Use as replacement frame for models 7150146 and 7150156 only
- Constructed of durable steel
- Must be used with locking band model 11274

Description

This sturdy metal frame (model 11273) is a replacement cradle for two of our Justrite steel dispensing cans (models 7150146 and 7150156 only). The frame allows the can to rotate into a stable horizontal position or setup vertically to tilt, making dispensing easy.

Note: This part number is for one replacement frame only. This does not include the locking band (model 11274) or a safety can.

Specifications

Model No	11273
UPC	697841004870
Color	Red
Material Specifications	Steel
Dimensions, Exterior	12.8 x 14.75 x 15.64
	32.5 x 37.5 x 39.7 cm