

# Add-On Flammable Safety Cabinets

## Product Images



## Short Description

---

- Flexible Installation: Can be mounted on top of existing safety cabinets, used as benchtop units, or wall-mounted for added flexibility
- Durable Construction: Built with sturdy 18-gauge double-wall, welded steel for maximum safety and compliance
- Secure Closure: Features a 3-point self-latching door system with a U-Loc™ handle, available in manual closing or self-closing options for enhanced safety
- Spill Protection: Includes a leakproof sump to capture dangerous spills and ensure safe storage
- Compliance and Durability: Meets OSHA & NFPA standards and is FM approved, with a durable lead-free epoxy/polyester powder-coat finish for chemical resistance
- Configurability: Available in four colors (yellow, red, gray, white) to meet specific storage needs
- Additional Storage: Adds up to 15 gallons of additional storage capacity

## Description

---

Eagle Add-On Flammable Cabinets offer a versatile solution for expanding storage space while safely segregating hazardous liquids. These cabinets can be mounted on top of existing 30- and 45-gallon Eagle safety cabinets or used independently as benchtop units, providing 15 gallons of additional storage capacity. They can also be mounted to a wall for added flexibility.

Constructed with sturdy 18-gauge double-wall, welded steel, these cabinets provide the same level of protection as taller safety cabinets, ensuring maximum safety and compliance. The 3-point self-latching doors close securely for maximum protection under fire conditions, and a leakproof sump captures dangerous spills. The U-Loc™ handle includes a double key set and accepts an optional padlock for increased security.

Eagle Add-On Safety Cabinets are designed to meet OSHA and NFPA standards and are FM-approved, ensuring compliance with industry regulations. To configure your cabinet, select from various color options, including yellow for flammables, red for combustibles, gray and white for waste storage. Choose between manual closing or self-closing doors to meet your specific storage needs.

Self-closing doors are a valuable feature in safety cabinets, offering enhanced safety and compliance. These doors automatically close and latch in the event of a fire or emergency, providing a critical barrier against the spread of fire and reducing the risk of ignition outside the cabinet. Self-closing doors minimize human error by ensuring that the cabinet is always properly closed, even when staff are distracted or forgetful. This feature is particularly important for storing flammable liquids, as it helps contain vapors and prevent accidental releases, thereby improving overall workplace safety. The self-closing mechanism uses fusible links that melt at 165°F (74°C), causing the doors to shut automatically under fire conditions.

These cabinets feature a durable lead-free epoxy/polyester powder-coat finish, providing excellent chemical resistance and minimizing the effects of corrosion and humidity. By configuring your Add-On Safety Cabinet, you ensure a safe and compliant storage solution that meets the unique requirements of your workplace.

Investing in Eagle Add-On Safety Cabinets not only enhances safety but also provides a flexible storage solution for hazardous materials, allowing you to optimize your workspace while adhering to industry standards.

## Specifications

---

Model No

ADD10

## Prop 65

---

**⚠️WARNING:** This product can expose you to chemicals including Lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

## Product Options

---

Color:

Red

Yellow

Gray

Beige

White

Door Style:

Self-close

Manual Close